Lab - 8



Sulphate

Gravimetric Method for Sulphate Determination:

Theory:

Nearly all sulphates in water can be precipitated as $BaSO_4$ on reaction with barium chloride ($BaCl_2$) under acidic conditions. Since this compound has very low solubility product, the precipitated $BaSO_4$ can be filtered and the residue weighed. The weight of the residue is proportional to the original sulphate concentration in solution.

Procedure:

- Dissolve 100 g of BaCl₂.2H₂O in 1000 mL of distilled water. Label the bottle, "Barium Chloride Solution for Sulphate Determination". (This solution will be available in the laboratory).
- Dissolve Na₂SO₄ in 1000 mL of water to make a sulphate solution containing 1000 mg sulphate/L. <u>Prepare four such bottles</u>. Label, "Sulphate Standard for Turbidimetric / Gravimetric Method". (This solution will be available in the laboratory).
- Add 50 mL water to 50 mL concentrated HCl. Label, "1+1 HCl solution". (*This solution will be available in the laboratory*).
- Take 100 mL aliquot each of sulphate standards of 0, 10, 25, 50 and 100 mg/L sulphate concentrations. Also prepare a 100 mL aliquot of sample.
- Adjust pH of each standard to 4.5-5 with the 1+1 HCL solution.
- ➢ While stirring gently, slowly add BaCl₂ solution until precipitation appears to be complete. Then add about 2 mL of excess BaCl₂ solution.
- Filter using GF/C filter, and dry at 103°C overnight. Cool in desiccator and weigh. Draw a calibration curve of the weight of residue versus sulphate concentration in solution.
- Use this calibration curve to determine the sulphate concentration in sample.

Turbidimetric Method for Sulphate Determination:

Theory:

Nearly all sulphates in water can be precipitated as BaSO₄ since this compound has very low solubility product. The precipitated BaSO₄ particles form a turbid suspension, turbidity of which is proportional to the original sulphate concentration in the solution.

Procedure:

- Buffer Solution: Dissolve 30 g MgCl₂.6H₂O, 5 g CH₃COONa.3H₂O, 1 g KNO₃, and 20 mL of acetic acid (99%)in 500 mL distilled water, and make up to 1000 mL. Label, "Buffer Solution for Sulphate Determination by Turbidimetric Method". (This solution will be available in the laboratory).
- Dissolve 1.479 g anhydrous Na₂SO₄ in 1000 mL of water to make a sulphate solution containing 1000 mg sulphate/L. Label, "Sulphate Standard for Turbidimetric/Gravimetric Method". (This solution will be available in the laboratory).
- Measure 100 mL of samples containing 0, 50, 100, 150, and 200 mg/L sulphate in 250 mL conical flasks. Add 20 mL of buffer solution and mix on a stirring apparatus. While stirring add a spoonful of BaCl₂ crystals and stir for a further 60 seconds.
- > Measure turbidity of the solutions for obtaining calibration curve.
- Measure sulphate concentration in sample using the above calibration curve.

8A Understanding About Sulphate in Water (4 x 10 Marks = 40 Marks)

- 1. What is the significance of a high sulphate concentration in water supplies and in wastewater disposal?
- 2. List four precautions that must be observed to ensure an accurate gravimetric determination of sulphate concentration.
- 3. What are two purposes for the conditioning reagent used in the turbidimetric determination of sulphate concentration?
- 4. From equilibrium considerations, calculate the relative proportions of sulphide in the H_2S , HS^- and S^{2-} forms at a) a pH of 6.0; b) a pH of 7.5; and c) a pH of 10.0.

8B Estimation of Sulphate in Water

A. Lab Report:

(4 x 10 Marks = 40 Marks)

Data: A Mass of unknown sample containing sulphate: 0.6039 g Mass of empty filter: 24.3849 g Mass of filter + BaSO4: 35.3639 g

Calculations (show work):

- 1. Calculate the mass of BaSO₄.
- 2. Calculate the mass of sulphate ion in the original solid.
- 3. Calculate the % mass of sulphate in the unknown sample.
- 4. The solid unknown was potassium sulphate. Calculate the percent error.
- **B.** The sulphate ion concentration in natural water can be determined by measuring turbidity that results when an excess BaCl₂ is added to a measured quantity of sample.

(1 x 4 Marks + 2 x 8 = 20 Marks)

A turbidimeter, the instrument used for the analysis, was calibrated with a series of standard NaSO₄ solutions. The following data were obtained in the calibration of sulphate concentration, C_x :

C _x , mg SO ₄ ²⁻ /L	Turbidimeter reading, R
0.00	0.06
5.00	1.48
10.00	2.28
15.00	3.98
20.00	4.61

Assume that a linear relationship exists between the instrument reading and the concentration.

- 1. Plot the data on a graph paper
- 2. Calculate using the least-squares method, slope and intercept, and write the equation of the line.
- Calculate the concentration of sulphate in a sample yielding a turbidimeter reading of 2.84. Find the standard deviation and the coefficient of variation.